

Dental instrument sterilization log

Check this template against your local regulations and your own policies before you use it.

Business _____ Location _____ Date _____

- Record every cycle, including daily test cycles and biological (spore) tests, as your infection control guidelines and the sterilizer's maker require.
- Answer Y or N for each check. Check the printout or display shows the cycle reached its time, temperature and pressure, and file the printout or data with this log.
- Check the chemical indicator on or in each pack has changed before the load is released, and write the cycle number on each pack.
- If a cycle or test fails, don't use the load. Reprocess it, record the failure, and follow your failure procedure before using the sterilizer again.

Date	Cycle no.	Sterilizer	Load contents	Parameters met?	Indicator passed?	Test result	Released by